## **MP RAPID 2019-NCOV IgG/IgM COMBO TEST CARD**

FOR THE QUALITATIVE ASSESSMENT OF IgG AND IgM ANTIBODIES TO  
2019 NOVEL CORONAVIRUS IN HUMAN SERUM, PLASMA, OR WHOLE BLOOD

**Catalog Number: 07RDCOV19**

***For In Vitro Diagnostic Use Only***

### INDEX

|          |                            |    |
|----------|----------------------------|----|
| English  | Instruction for use        | 2  |
| Deutsch  | Gebrauchsanweisung         | 6  |
| Français | Instructions d'utilisation | 10 |
|          | Symbols                    | 14 |

**INTENDED USE**

MP Rapid 2019-nCoV IgG/IgM Combo Test Card is an immunochromatography based one-step in vitro test. It is designed for the rapid qualitative determination of IgG and IgM antibodies to 2019 novel coronavirus (2019-nCoV, SARS-CoV-2) in human serum, plasma, or whole blood. MP Rapid 2019-nCoV IgG/IgM Combo Test Card is a supplementary detection for COVID-19 suspected infected patients besides nucleic acid testing, which could greatly raise the accuracy of the detection for COVID-19.

**SUMMARY**

Corona Virus Disease 2019 (COVID-19) is an acute infectious disease caused by 2019 novel coronavirus (SARS-CoV-2). The incubation period of the disease is 1-14 days, mostly 3-7 days. The incubation period is infectious. Asymptomatic infections may also be the source of infection. Respiratory droplets and contact are the main routes of transmission. The initial symptoms of patients are fever, fatigue and cough and patients can also gradually develop dyspnea and other serious manifestations. Most patients have a good prognosis, however some severe cases may develop acute respiratory distress syndrome, septic shock or even death. At present, there is no specific treatment for the disease.

IgM antibodies can be detected soon after the incubation period and remain present for a short time. IgM positive blood samples can be an indicator of an acute infection. IgG antibodies appear a few days into the incubation period and remain present for a long time. IgG positive blood samples can be an indicator of an infected or previously infected patient.

**PRINCIPLE**

MP Rapid 2019-nCoV IgG/IgM Combo Test Card utilizes the principle of immuno-chromatography. Mouse anti-human IgM and human IgG antibodies are immobilized on the nitrocellulose membrane respectively, as two individual test lines (IgM line and IgG line) in the test window of the test device. The IgM line in the test window is closer to the sample well, and is followed by the IgG line. As the test sample flows through the membrane within the test device, the colored-2019-nCoV antigen-colloidal gold conjugate complexes with specific antibodies (IgM and/or IgG) of 2019 novel coronavirus, if present in the sample. This complex moves further on the membrane to the test region where it is captured by the anti-human IgM and/or human IgG antibodies coated on the membrane leading to formation of a colored band, which indicates a positive test result. Absence of this colored band in the test window indicates a negative test result. A built-in control line will always appear in the test window when the test has performed properly, regardless of the presence or absence of anti-2019 novel coronavirus antibodies in the specimen.

**MATERIALS PROVIDED**

1. MP Rapid 2019-nCoV IgG/IgM Combo Test Card
2. Sample buffer
3. 2 µL capillary pipet
4. Instructions for Use

**MATERIALS REQUIRED BUT NOT SUPPLIED**

Clock or timer, safety lancets, alcohol prep-pad, specimen collection container, centrifuge, biohazard waste container, disposable gloves, disinfectant.

**STORAGE**

1. Store the test device at 2 to 30°C in the original sealed pouch. Do Not Freeze.
2. The expiration date indicated on the pouch was established under these storage conditions.
3. The test device should remain in its original sealed pouch until ready for use. After opening, the test device should be used immediately. Do not reuse the device.

**PRECAUTIONS**

1. For **professional *in vitro*** diagnostic use only.
2. Do not use the product beyond the expiration date.
3. Do not use the product if the pouch is damaged or the seal is broken.
4. Handle all specimens as potentially infectious.
5. Follow standard laboratory procedures and biosafety guidelines for handling and disposal of potentially infectious material. When the assay procedure is completed, dispose specimens after autoclaving at 121° C for at least 20 min or treating with 0.5% Sodium Hypochlorite for 1-2 hours.




**SPECIMEN COLLECTION AND PREPARATION**

1. The serum, plasma or whole blood specimen should be collected under standard laboratory conditions.
2. Heat inactivation of specimens, which may cause hemolysis and protein denaturation, should be avoided.
3. Patient samples have the best performance when tested immediately after collection. If specimens are to be stored, the red blood cells should be removed to avoid hemolysis. If the assay is not performed immediately, serum or plasma specimens may be refrigerated at 2-8°C up to 3 days. For long-term storage, freeze the specimen at -20°C for 3 months or -70°C for longer periods. Allow sample to reach room temperature before proceeding.
4. Sodium azide can be added as a preservative up to 0.1% without affecting the test results.

**QUALITY CONTROL**

1. The control band is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagents are reactive.
2. Good Laboratory Practice recommends the daily use of control materials to validate the reliability of the device. Control materials which are not provided with this test kit are commercially available.

**PROCEDURE**

|  |  |   |
|--|--|---|
| <b>1</b>   | Bring the kit components to room temperature before testing.                             |   |
| <b>2</b>   | Open the pouch and remove the Card. Once opened, the test card must be used immediately. |   |
| <b>3</b>   | Label the test card with patient identity.   |   |
| <b>4</b>   |  |   |
|      |        | Withdraw the blood specimen with the capillary pipet provided. Gently squeeze out the extra specimen to leave 2µL in the pipet as marked with the scale line. Apply 2µL of blood specimen to the “S1” area as marked. |
| <b>5</b>   | <b>6</b>   |   |
|     |      |   |
| Add 2 drops of sample buffer (approximately 80-100µL) to the well which is marked “S”. |  | Read the result at 15 minutes. A strong positive sample may show a result earlier.<br><b>Note: Results after 20 minutes may not be accurate</b>   |

**INTERPRETATION OF RESULTS**

| <b>POSITIVE</b>  |  |   |
|--|--|---|
|                       |   |    |
| <b>Both IgG/IgM Positive</b>   | <b>IgM Positive<br/>IgG Negative</b>   | <b>IgM Negative<br/>IgG Positive</b>  |
| Control line and both test lines appear. It indicates the possibility of an acute secondary infection. | Both control line and the second test line (the lower test line which is closer to the sample well) appear. It indicates the possibility of a primary infection. | Both control line and the first test line (the higher test line) appear. It indicates the possibility of a secondary infection or a past infection. |
| <b>NEGATIVE</b>  |  |   |
|                      | Only control line appears.   |   |
| <b>INVALID</b>   |  |   |
|                     | The test result is invalid if a colored band does not form in the control region. The sample must be re-tested, using a new test device.                         |   |

**PERFORMANCE CHARACTERISTICS:****Accuracy**

A total of 74 specimens from confirmed patients were tested, the results showed that 65 specimens were IgM positive and / or IgG positive, and the clinical sensitivity was 87.8%. A total of 305 specimens from healthy persons were tested, the results showed that 302 specimens were both IgM and IgG negative, 1 specimen was IgM positive, 2 specimens were IgG positive, and the clinical specificity was 99.0%. The accuracy was 96.8%.

**Assay Specificity****1. Other infectious diseases**

MP Rapid 2019-nCoV IgG/IgM Combo Test has tested samples that were infected by the following diseases: Influenza A Virus, Influenza B Virus, Adenovirus, Rotavirus and Mycoplasma Pneumoniae. All the samples showed no effect on the specificity of the assay.

**2. Blood compounds**

MP Rapid 2019-nCoV IgG/IgM Combo Test has tested samples with high Rheumatoid Factor (RF), Bilirubin, Triglyceride and Hemoglobin. The results showed that these compounds had no effect on the specificity of the assay up to the listed concentration.

|                   |            |
|-------------------|------------|
| Rheumatoid Factor | 80 IU/ml   |
| Bilirubin         | 342 µmol/L |
| Triglyceride      | 37 mmol/L  |
| Hemoglobin        | 10 mg/mL   |

**3. Common drugs**

MP Rapid 2019-nCoV IgG/IgM Combo Test has tested samples with common drugs. The results showed that these drugs had no effect on the specificity of the assay.

Histamine Hydrochloride, Interferon- $\alpha$ , Zanamivir, Ribavirin, Oseltamivir, Peramivir, Lopinavir, Ritonavir, Arbidol, Levofloxacin, Azithromycin, Ceftriaxone, Meropenem, Tobramycin.

**LIMITATIONS**

1. The test is limited to the qualitative detection of anti-2019-nCoV antibody levels in serum, plasma, or whole blood specimens. The exact concentration of anti-2019-nCoV antibody cannot be determined by this assay.
2. Although the test is very accurate in detecting anti-2019-nCoV antibodies, a low incidence of false results can occur. Other clinically available tests are required if questionable results are obtained.
3. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
4. In the early stage of infection, if IgM and IgG antibodies are not produced or the titer is very low, false negative results will be observed. It is suggested that patients should collect samples again after 7-14 days, and test along with the previously collected samples to confirm whether there is serological positive transfer or significant increase in titer. In the later stage of infection, IgM titer will decrease or even be negative, while IgG will continue to increase.

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Manufactured by/ Hergestellt durch/ Fabricant: MP Biomedicals Germany GmbH